

# Conquer O2

Professional encapsulating oxidising multi-task spotter



Conquer O2 has been designed to build on the outstanding ability of free oxygen to remove the toughest of stains from fibres & fabrics, provide powerful sanitisation and eliminate mal-odours. Conquer O2 has teamed up the awesome oxidising power of oxygen with

wetting agents and encapsulating polymers to deliver a carpet and fabric stain remover which will remove those impossible stains which other spotters leave behind.

Thanks to the latest in polymer and surfactant technology, Conquer O2 delivers reliable, dependable results every time.

## BENEFITS

- **Stain Removal** – Lifts stubborn stains caused by foods, drinks, cosmetics, organics, beverages, pets.
- **Blood, Urine & Body fluids** – Conquer O2 provides powerful stain removal, sanitizing and deodorising action on blood, urine, body fluids and other organics.
- **Versatile** – Because Conquer dries to a brittle residue, it makes an ideal product to clean soil from garments, drapes and soft furnishings or partitions.
- **Slows resoiling** – Advanced polymer prevents rapid resoiling.
- **Non-rinsing** – Soil residues dry to a brittle residue for vacuuming away.

## DIRECTIONS for USE

Always pretest in an inconspicuous location for fibre and dye sensitivity.

1. Ready to use – no dilution needed.
2. Sponge or scrape up any surface spot material.
3. Spray or blot solution onto the spot. Do not overwet. Agitate gently.
4. Allow 15 minutes dwell time for stubborn stains.
5. Blot away with a clean white towel. Do not rub. Stubborn spots may need light agitation. Work from the outside to the centre of large spots
6. Repeat as necessary.

Attend to spots and spills as soon as possible. Fresh spots are significantly easier to remove. If left they can become permanent stains.

## PRODUCT SPECIFICATIONS

**Product code:** AP478

**pH:** 5.5 RTU

**Dilution:** 1:1 for delicate fibres

**Fragrance:** None

**Safety:** Non-hazardous

**Colour:** Clear liquid

**Packaging:** 500ml, 5Lt & 20Lt

